

ISDH Radiologic Health Section  
ISDH THERAPY ACCELERATOR  
INITIAL COMMISSIONING SUMMARY

1. Facility Registration Number \_\_\_\_\_ 2. Physicist's Signature \_\_\_\_\_  
3. Inspection Date \_\_\_\_\_ (code)

A. Machine Number	B. Machine Design (Use Codes)	C. Location	D. Manufacturer (Use Codes)	E. Date of Manufacturer	F. Model and Serial Numbers

**RADIATION SAFETY SURVEY:**

Qualified Physicist who performed shielding calculations: Code \_\_\_\_\_ Shielding Document Submitted: ☐ Yes ☐ No  
Workload Assumptions: \_\_\_\_\_

Measurement Set-up: \_\_\_\_\_

Radiation Safety Instrumentation / Calibration Date: \_\_\_\_\_

Radiation Survey Results	mRem	mRem	mRem	mRem	
Location <sup>1</sup>	0°	90°	180°	270°	mRem / week
Console / Control Area					
Primary Barrier (90° Gantry Angle)					
Primary Barrier (270° Gantry Angle)					
Roof					
Door					
Secondary Barrier					
Secondary Barrier					
Secondary Barrier					
Secondary Barrier					
HVAC Ductwork					
Accessible Conduits					

<sup>1</sup> Submit a facility layout indicating where the measurements were taken.

**DOSIMETRY SYSTEM AT FACILITY:**

Cylindrical Chamber: Manufacturer/Model \_\_\_\_\_ ADCL Calibration Date: \_\_\_\_\_  
 Parallel Plate Chamber: Manufacturer/Model \_\_\_\_\_ N<sub>gas</sub> Derivation Date: \_\_\_\_\_  
 Electrometer: Manufacturer/Model \_\_\_\_\_ ADCL Calibration Date: \_\_\_\_\_  
 Date of Aneroid Barometer Intercomparison: \_\_\_\_\_ Date Initial Calibration Completed: \_\_\_\_\_  
 Calibration Protocol: \_\_\_\_\_

**QUALITY MANAGEMENT PROGRAM:**

Submit the following departmental quality assurance documents:

- Treatment planning computer and dose calculational algorithm QA procedures
- Patient chart review policy and procedures
- Weekly output constancy check policy and procedure
- Monthly output spot check procedure

**FACILITY REQUIREMENTS:**

Warning Light at Entrance to Accelerator Room Operational Sat ☐ Unsat ☐

**ACCELERATOR REQUIREMENTS:**

Absorbed Dose due to Leakage Radiation in the Patient Plane:	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Transmission Through Collimators:	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Removable Wedges Clearly Identified:	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
X-ray Contamination of Electrons Beams:	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/> Not Applicable <input type="checkbox"/>
Absorbed Dose at the Surface:	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Dose Monitoring System:		
Incorporated into two independent dose monitoring systems	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Monitor units displayed in the event of a power failure	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Asymmetry of beam monitored and interlocked	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Secondary dose monitoring system able to terminate irradiation	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Timer terminates irradiation if dose monitoring systems fails	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/>
Selection of Radiation Mode	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/> Not Applicable <input type="checkbox"/>
Selection of Energy	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/> Not Applicable <input type="checkbox"/>
Selection of Fixed or Rotational Mode	Sat <input type="checkbox"/>	Unsat <input type="checkbox"/> Not Applicable <input type="checkbox"/>

Mechanical Parameters		Tolerance <sup>2</sup>
Gantry Rotation Isocenter	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Collimator Rotation Isocenter	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Table (Couch) Rotation Isocenter	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Field Size Readouts	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Gantry Angle Readouts	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Collimator Angle Indicator	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
ODI Accuracy & Linearity	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Light/Radiation Field Congruency	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Laser Alignment System Accuracy	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	

  

Emergency Off Switches	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>
Door Interlock	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>
Dosimetry Interlocks <sup>3</sup>	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>
Safety Interlocks <sup>4</sup>	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>
Audio Communications System	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>
Visual Monitor System	Operational: Yes <input type="checkbox"/> No <input type="checkbox"/>

Photon Beam Parameters				Tolerance <sup>2</sup>
Nominal Energy (MV)				-----
Ionization Ratio				-----
Flatness	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Symmetry	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Field Size Factors	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
PDD / TPR Values	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Output Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed
Output / MU Reproducibility	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Monitor Chamber Linearity	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Wedge / Compensator Transmission Factors	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Tray / Custom Block Transmission Factors	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Open Field Isodose Lines (central axis)	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Wedge Isodose Lines (central axis)	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Constancy Check Device Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed

<sup>2</sup>The qualified radiation oncology physicist shall specify tolerance values based upon accelerator manufacture specifications and the most recent published standards. Corrective action is required for measured data that exceeds the stated tolerance.

<sup>3</sup>Energy (bending magnet current), flatness, symmetry, temperature and compensation, et cetera.

<sup>4</sup>Target slide or scattering foil, dose chamber slide, dose rate, et cetera.

Electron Beam Parameters				Tolerance <sup>2</sup>
Nominal Energy (MeV)				-----
Mean Incident Energy				-----
Flatness: (Cone Size _____)	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Symmetry	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Output Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed
Output / MU Reproducibility	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Constancy Check Device Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed

Electron Beam Parameters				Tolerance <sup>2</sup>
Nominal Energy (MeV)				-----
Mean Incident Energy				-----
Flatness: (Cone Size _____)	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Symmetry	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Output Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed
Output / MU Reproducibility	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Constancy Check Device Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed

Electron Beam Parameters				Tolerance <sup>2</sup>
Nominal Energy (MeV)				-----
Mean Incident Energy				-----
Flatness: (Cone Size _____)	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Symmetry	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Output Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed
Output / MU Reproducibility	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	Sat <input type="checkbox"/> Unsat <input type="checkbox"/>	
Constancy Check Device Calibration	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>	Performed